



## Fronts and extension tips

# Fronts metal

nickelfree

## Model „Nikolai“

Size	Measurement		Article Number					
	B	C	Gold	Silver	Black	Grey	Red	Blue
55 / 20	55 mm	20 mm	77733170	77733171	77733172	77734784	77734785	77734783
57 / 20	57 mm	20 mm	77733167	77733168	77733169	77734787	77734788	77734786
59 / 20	59 mm	20 mm	77733381	77733382	77733383			



Colour: Gold



Colour: Silver



Colour: Black



Colour: Grey



Colour: Red



Colour: Blue

## Model „Sportivo“

Size	Measurement		Article Number	
	B	C	Silver	Black
59 / 14	59 mm	14 mm	77733179	77733352
63 / 14	63 mm	14 mm	77733386	77733387



Colour: Silver



Colour: Black

# Fronts metal

## Model „Josephine“

Size	Measurement		Article Number					
	B	C	Gold	Silver	Black	Grey	Red	Blue
51 / 18	51 mm	18 mm	77733375	77733376	77733377	77734790	77734791	77734789
53 / 18	53 mm	18 mm	77733176	77733177	77733178	77734793	77734794	77734792
55 / 18	55 mm	18 mm	77733173	77733174	77733175			



Colour: Gold



Colour: Silver



Colour: Black



Colour: Grey




Colour: Red



Colour: Blue

## Model „Lars“

Size	Measurement		Article Number
	B	C	Silver
56 / 16	56 mm	16 mm	77734602



Colour: Silver

# Fronts plastic



## Model „Thomas“

Size	Measurement		Article Number
	B	C	
53 / 20	53 mm	20 mm	77733949
59 / 20	59 mm	20 mm	77734428



## Model „Erik“

Size	Measurement		Article Number
	B	C	
53 / 20	53 mm	20 mm	77733950
59 / 20	59 mm	20 mm	77734427



## Model „Celine“

Size	Measurement		Article Number
	B	C	
51 / 20	51 mm	20 mm	77734429




## Model „Frances“

Size	Measurement		Article Number
	B	C	
51 / 20	51 mm	20 mm	77734430
53 / 20	53 mm	20 mm	77734885



## Model „Florence“

Size	Measurement		Article Number
	B	C	
51 / 20	51 mm	20 mm	77734426
53 / 20	53 mm	20 mm	77733933



Nickelfree

# Extension tips metal

for plastic fronts: Florence, Thomas, Erik

Size	AN-Evo1 Total length appr.		contact star evo1 Total length appr.		Article Number Silver
	mm	inch	mm	inch	
2	82,00	3,22	77,00	3,02	77734431
3	85,50	3,36	80,50	3,16	77734110
4	89,00	3,50	84,00	3,30	77733940
5	92,50	3,64	87,50 <td 3,44	77733941	
6	96,00	3,77	91,00	3,57	77733942



Colour: Silver




Fig. 1  
To determine the exactly length of extension tips, please use measure tool MA1. (See accessories)




Fig. 2  
Mount the cover to extension tip with 2 pieces screw No. 2. Then mount the complete extension tip on the frame with screw No. 1.




Fig. 3  
After snap-on the bone conducted spectacles you will get the approximately length of the temple.



# Extension tips

www.bhm-tech.at

made of metal: Nikolai, Josephine, Sportivo, Lars, Celine, Frances

Size	Total length approximately				Article Number					
	AN-Evo1		contact star evo1		Gold	Silver	Black	Grey	Red	Blue
	mm	inch	mm	inch						
1	78,50	3,08	73,50	2,89	77733367	77733369	77733368			
2	82,00	3,22	77,00	3,02	77733356	77733357	77733358			
3	85,50	3,36	80,50	3,16	77733184	77733187	77733190	77734796	77734797	77734795
4	89,00	3,50	84,00	3,30	77733185	77733188	77733191	77734799	77734800	77734798
5	92,50	3,64	87,50	3,44	77733186	77733189	77733192	77734802	77734803	77734801
6	96,00	3,77	91,00	3,57	77733396	77733397	77733398			
7	99,50	3,91	94,50	3,71	77733399	77733400	77733401			
8	103,00	4,05	98,00	3,85	77733402	77733403	77733404			



Colour: Gold



- Colour: Silver
- Colour: Black
- Colour: Red
- Colour: Blue
- Colour: Grey
- Colour: Gold

## For your individual mounting

Hearing system	Size	Total length approximately	Article Number	
			Brown	Black
contact star evo1	XL	up to 108 mm	77734433	77734434
AN-Evo1	XL	up to 118 mm	77733196	77733197

Description	Article Number
coupler for temples L+R	77734876



Colour: Brown



Colour: Black



# Accessories

## Presentation case

Description	Measurement			Article Number
	Length	Width	Height	
Presentation case made of plastic	340 mm	300 mm	70 mm	77733546



## Frame heater



In order to adapt (bend) the hearing aid spectacles, it must first be heated in the bending area designed for the purpose. Further information is available in our folder „Tips & Tricks“

Please don't hesitate to contact us if you require a frame heater!

## Measure tool

Description	Article Number
Measure tool MA1	77733215





BHM-Tech Produktionsgesellschaft mbH  
Grafenschachen 242  
7423 Grafenschachen, Austria  
[www.bhm-tech.at](http://www.bhm-tech.at)



According to EU directive 93/42/EEC

**0297** Quality management system according to DIN EN ISO 13485

Perfection made in Austria

