

Bone conduction hearing system

	contact mini	bone conduction implant device
Purpose / Intended use 	Provides a comfortable solution to aid in bone conduction hearing. Designed for incorporated use in a wide range of headgears. It is worn externally with no surgery required.	Helps in bone conduction hearing. The device can be worn with a soft band. Surgery will be required in the long run to get the maximum benefit.
Technical details	contact mini	bone conduction implant device
feedback cancellation	✓	✓
noise reduction	✓	✓
multi channel	✓	✓
tele coil integrated in device	✓	✗ (external part)
Wearing possibilities	contact mini	bone conduction implant device
softbands 	✓ (3 colours available)	✓
headbands 	✓ (different colours and sizes available)	✗
sports caps 	✓ (2 colours available)	✗ (due to feedback sensitivity)
hairbands 	✓ (different sizes available)	✗
invisible device (decorative elements for covering) 	✓ (different motifs available)	✗ (due to feedback sensitivity)
individual headpiece solutions (integration in sport caps, helmets, ... on demand)	✓	✗
additional head pieces can be worn (e.g. usage of sun hat, ...)	✓	✗ (due to feedback sensitivity)
glasses can be worn at the same time	✓	✗ (due to feedback sensitivity)
flexible wearing options (headpieces can be exchanged)	✓	✗