



## Bone conduction hearing system

	contact mini	bone conduction implant device
<b>Purpose / Intended use</b> 	Provides a comfortable solution to aid in bone conduction hearing. Designed for incorporated use in a wide range of headgears. <b>It is worn externally with no surgery required.</b>	Helps in bone conduction hearing. The device can be worn with a soft band. <b>Surgery will be required in the long run to get the maximum benefit.</b>
<b>Technical details</b>	<b>contact mini</b>	<b>bone conduction implant device</b>
feedback cancellation	✓	✓
noise reduction	✓	✓
multi channel	✓	✓
tele coil integrated in device	✓	✗ (external part)
<b>Wearing possibilities</b>	<b>contact mini</b>	<b>bone conduction implant device</b>
<b>softbands</b> 	✓ (3 colours available)	✓
<b>headbands</b> 	✓ (different colours and sizes available)	✗
<b>sports caps</b> 	✓ (2 colours available)	✗ (due to feedback sensitivity)
<b>hairbands</b> 	✓ (different sizes available)	✗
<b>invisible device</b> (decorative elements for covering) 	✓ (different motifs available)	✗ (due to feedback sensitivity)
<b>individual headpiece solutions</b> (integration in sport caps, helmets, ... on demand)	✓	✗
<b>additional head pieces can be worn</b> (e.g. usage of sun hat, ...)	✓	✗ (due to feedback sensitivity)
<b>glasses can be worn at the same time</b>	✓	✗ (due to feedback sensitivity)
<b>flexible wearing options</b> (headpieces can be exchanged)	✓	✗